

TG extraction from seeds

Day 1

Lipid extraction by Folch's method; solvent evaporation; weighing; storage at -30°C.

Day 2

Separation of lipids by TLC; recovery of lipid spots; antioxidant addition; HCl/methanol treatment for GLC.

Days 3 and 4

Hexane extraction; GLC analysis.

Lipid extraction from Arabidopsis seeds

Sampling Table

Sample name (date):

Sample No.	Example	No. 1	No. 2	No. 3
Fresh weight (mg)	10			
H ₂ O from material (ml) — (W)	Ignore	Ignore	Ignore	Ignore
2-Propanol (ml)	1	1	1	1
<u>Homogenate (ml)</u> Methanol	0	0	0	0
Chloroform (1)	2	2	2	2
<u>Pellet washing (ml)</u> Methanol	1	1	1	1
Chloroform (2)	2	2	2	2
<u>Combined supernatant</u> (A + B in the text)	6	6	6	6
0.9% KCl	1.2	1.2	1.2	1.2
Empty vessel (g)				
Total lipid extracts (mg)				

Day 1

Place 1 ml 2 propanol in a 10 ml screw-capped glass tube and keep the tube in a 90 °C water bath.

↓

Add weighed seeds (10 mg) to the tube and keep for another 5 min.

↓

Chill the tube on ice.

↓

Add 2 ml chloroform, and homogenize the seeds, using a Polytron (max. speed).

↓ 3000 rpm x 15 min (r.t.)

Recover the supernatant (A) into a fresh 10-ml screw-capped glass tube by decantation (if the pellet is soft, remove by pipetting).

↓

Resuspend the pellet in 1 ml methanol and 2 ml chloroform (as listed in the **Sampling Table**) and mix well using a vortex mixer.

↓ 3000 rpm x 15 min (r.t.)

Recover the supernatant (B) and combine with the supernatant A.

↓

Add 1.2 ml (0.2 vol) of 0.9% KCl and mix well by vigorous shaking.

↓ 3000 rpm x 15 min (r.t.)

Recover the lower (chloroform) layer*¹ into a 20 ml weighed pear-shaped flask, and evaporate the solvent on a rotary evaporator*².

*¹Special attention should be paid not to recover any material floating in the intermediate phase.

*²If you see water droplets on the bottom of flask, add a few ml of ethanol and evaporate again.

↓

Dry up the residual solvent in a vacuum desiccator and determine the lipid dry weight.

↓

Add chloroform to make a 10 mg/ml lipid solution while placing the flask on ice to prevent solvent evaporation.

↓


Transfer the lipid solution into a 1-ml screw-capped sample tube*³.

*³Because the cap of the tube is so small, special attention should be paid to ensure that the cap is screwed tightly on the tube. Each cap should then be sealed with parafilm and stored at -30 °C. If you store the sample for more than one week, the antioxidant BHT should be added at a final concentration of 10 µg/ml.

Day 2

Separation of lipids by silica gel thin-layer chromatography

Using a 100- μ l glass micro syringe, streak 100 μ l of the lipid extract (1 mg lipid equivalent) on a 1-cm line marked 2 cm from the bottom of TLC plate (20 x 20 cm, Silica gel 60, Merck 1.05721.) as shown below.



1st dimension :

Hexane/diethyl ether/acetic acid = 80:30:1 (by volume)



Air-dry the plate (you may use warm air from a hair dryer).



Spray 0.01% Purimurin/80% acetone solution over the plate.

Dry the plate using warm air from a hair dryer.



Mark lipid spots under long-wave UV light (360 nm), using a soft pencil (HB).
The fluorescence image must be taken on a digital camera.



Scrape off each spot on medicine paper, and carefully transfer all silica gel powder into a 15-ml screw-capped glass tube^{*4}.

^{*4}Special attention should be paid not to use any tubes with cracks on the mouth or the screw mound.



Add 10 µl each of a 15:0 standard/toluene solution, using a fresh glass capillary pipette.

↓

Add 10 µl of a BHT/toluene solution, using a fresh glass capillary pipette.
(If 14:0 is expected in the lipid samples, BHT must not be added, because 14:0-Me and BHT show the same retention time on a gas chromatograph.)

↓

Add 3 ml each of HCL/methanol solution. Screw on the tube and mix vigorously with a vortex mixer.

↓ 80 °C for 3 h on an aluminum heating block.

Day 3

Chill all tubes on ice.



Extract FA methyl esters with 3 ml distilled hexane.



Evaporate hexane, using a vacuum centrifugation evaporator. Dissolve dried lipid residues in ~5 drops of hexane from a Pasteur pipette and subject to gas chromatograph analyses.

(The hexane solution should be transferred to a 1-ml sample tube for storage or a glass tube for an automatic sample injector.)

Correspondence: Ikuo Nishida, nishida@molbiol.saitama-u.ac.jp