AtMap1 DNA microarray experiment

- 1. > 100 seedlings are collected and lyophilized.
- 2. Genomic DNA is extracted from the lyophilized materials with a DNeasy Plant Mini Kit (Qiagen, Valencia, CA).
- 3. Double-stranded DNA is quantified with a Quant-iT dsDNA Broad-Range Assay Kit (Invitrogen, Carlsbad, CA).
- 4. The genomic DNA (1 μg) is used for labeling reaction with a Genomic DNA Labeling Kit (Agilent technology) according to direct labeling protocol, except following points.
- The fragmentation of genomic DNA is performed with two sets of restriction enzymes instead of the sonication.
- A half of the genomic DNA (0.5 μg) is digested with RsaI/AluI and the rest (0.5 μg) is digested with DraI.
- 5. Hybridization is performed with an Oligo aCGH/ChIP-on-Chip Hybridization Kit (Agilent Technology) according to the standard protocol, except that **Cot-1 DNA** is not added to the hybridization solution.
- 6. After hybridization, the microarray is washed in Oligo aCGH/ChIP-on-Chip Wash Buffer 1 and 2 (Agilent Technology). Immediately after washing the microarray slide scanned by a DNA microarray scanner (Agilent Technologies). The signals are quantified and normalized with Feature Extraction Software ver.
- 9.1 (Agilent Technologies) using the software's default parameters.

- 7. The design file of AtMap1 for the Feature Extraction is available from http://podb.nibb.ac.jp/Organellome/FunctionalPDF/MAP1/ExFullGEML.xml.
- 8. The annotation file for AtMap1 is available from http://podb.nibb.ac.jp/Organellome/FunctionalPDF/MAP1/ExFullGEML.xml. The file contains information about the features of each probe (gene, intergenic region and transposable element), the AGI code of the nearest gene and the annotation about the gene (TAIR8).
- 9. The .BED design file of AtMap1 for the UCSC Genome Browser is available from http://podb.nibb.ac.jp/Organellome/FunctionalPDF/MAP1/015638_BED.bed.

Reference

Nagano, A.J., Fukazawa, M., Hayashi, M., Ikeuchi, M., Tsukaya, H., Nishimura, M., and Hara-Nishimura, I. (2008). AtMap1: a DNA microarray for genomic deletion mapping in *Arabidopsis thaliana*. Plant J. 47, in press. (doi: 10.1111/j.1365-313X.2008.03656.x)