

Preparation of crude membrane fractions from rice tissues and organs

1) Solutions

Extraction buffer

50 mM KH ₂ PO ₄ -KOH, pH 7.5 at 4°C	Stock soln. 500 mM (1/10 dilution)
2 mM EGTA	Stock soln. 50 mM (1/25 dilution)
1 mM EDTA	Stock soln. 100 mM (1/100 dilution)
2 mM DTT	Conc. soln. 100 mM (1/50 dilution)
20 µM leupeptin	Conc. soln. 2 mM (1/100 dilution)
500 µM AEBSF	Conc. soln. 10 mM (1/20 dilution)
0.25 M D(-)-Mannitol	0.04554 g/mL
0.1 % (w/v) BSA	1 mg/mL

*AEBSF: 4-(2-aminoethyl)-benzenesulfonyl fluoride, Sigma A8456, AEBSF is irreversible serine protease inhibitor with high water solubility.

D1 buffer

10 mM KH ₂ PO ₄ -KOH, pH 7.8 at 4°C	Stock soln. 500 mM (1/50 dilution)
1 mM DTT	Conc. soln. 100 mM (1/100 dilution)
0.25 M D-Sorbitol	0.04554 g/mL

Suspension buffer

25 mM Hepes-BTP, pH 7.2 at 4°C	Stock soln. 250 mM (1/10 dilution)
* BTP: BisTris Propane	
1 mM DTT	Conc. soln. 100 mM (1/100 dilution)
500 µM AEBSF	Conc. soln. 10 mM (1/20 dilution)
0.25 M D-Sorbitol	0.04554 g/mL

2) Preparation of crude membrane fraction

On ice

One to 2 g f.w. of tissues or organs harvested.

↓

Add 1 to 2 ml of ice-cooled extraction buffer per 1 g f.w. of samples.

Homogenize by using a mortar and pestle.

↓

Centrifuge at 10,000 xg for 20 min at 4°C.

Collect the supernatant.

↓

Divide the supernatant into two Open Top Thick-walled polycarbonate tubes (Beckman; capacity: 2 ml) and fill up with the extraction buffer.

Set tubes on TLS-55 rotor (Beckman; Swing-type) pre-cooled at 4°C after balanced.

Set the rotor on Optima TLX ultracentrifuge (Beckman).

Centrifuge at 156,000 xg (50,000 rpm) for 30 min at 4°C.

↓

Transfer the supernatant as the soluble protein fraction into the new tube.

Recover the precipitation as the crude membrane fraction.

↓

Add 1 mL of the D1 buffer to the precipitation.

Re-suspend the precipitation by pipetting.

Fill up with the D1 buffer.

↓

Re-centrifuge at 156,000 xg (50,000 rpm) for 30 min at 4°C.

↓

Recover the precipitation as the crude membrane fraction.

↓

Added 10 to 40 µL of the suspension buffer to the prep.

Re-suspend the precipitation by pipetting.

↓

Quantify protein amounts in soluble and crude membrane fractions by the Bradford method using BSA as a standard.

↓

Divide into aliquots.

Freeze quickly with liquid nitrogen and stored at -80°C until used.