

Fixation


Cotyledons and Leaves

Fixative*¹

4% paraformaldehyde
1% glutaraldehyde
0.05M cacodylate buff. (pH7.4)

Washing buff.

0.05M cacodylate buff. (pH7.4)



*² Cut slightly above the base of the petioles

*¹ In case of etiolated cotyledons, use the fixative for dry seeds


Cut off the petioles (slightly above the base)*² and soak the samples in fixative*³

↓
Vacuum infiltrate in fixative
(about -99kPa 10min → slowly return to normal pressure, 3 times)

↓
Cut the samples into 1-2mm square sections

↓
Fix for 2hours in fixative (room temp.)

↓
Wash twice with washing buff.



*³ When the samples float...

Dry seeds

Fixative

10% dimethylsulfoxide
4% paraformaldehyde
1% glutaraldehyde
0.05M cacodylate buff. (pH7.4)

Washing buff.

0.05M cacodylate buff. (pH7.4)

Cut the seeds in half and soak them in fixative



Vacuum infiltrate

(about -99kPa 10min → slowly return to normal pressure, 3 times)



Fix for 2hours in fixative (room temp.)



Wash twice with washing buff.

Roots

Fixative

4% paraformaldehyde

1% glutaraldehyde

0.05M cacodylate buff. (pH7.4)

Washing buff.

0.05M cacodylate buff. (pH7.4)

Cut the samples into 1-2mm long sections and soak them in fixative



Vacuum infiltrate

(about -99kPa 10min → slowly return to normal pressure, 3 times)



Fix for 2hours in fixative (room temp.)



Wash twice with washing buff.

Embedding

Dehydration and infiltration

Replace the buffer with 50% (v/v) DMFA (Dimethylformamide)

50% DMFA 4°C 15min x 2

80% DMFA -20°C 15min x 2

90% DMFA -20°C 15min x 2

100% DMFA -20°C 20min x 3

LR White : DMFA = 1:1 -20°C 30min

LR White : DMFA = 2:1 -20°C 30min

LR White -20°C 30min

LR White -20°C overnight

Polymerize

Fill the "Beem capsules" with LR White (accelerator added)

Put the samples in the capsules

Close the capsules without leaving any air inside^{*4}

UV irradiation -20°C 24hrs

UV irradiation 4°C 1hr

UV irradiation room temp. 1-3hrs

^{*4} Add the resin until it overflows before closing the capsules

Sectioning

Trim the polymerized blocks and make ultrathin sections using ultramicrotome

Mount the sections on nickel grids

Immunostaining

Phosphate buffer saline (PBS)

0.01M phosphate buff.–0.85% NaCl

Blocking solution

1% BSA–0.1% NaN₃ in PBS

(All antibodies are diluted with blocking solution)

Blocking

Room temp. 30min in blocking solution

Wash grids with PBS

1st antibody reaction

4°C overnight in 1st antibody solution (in moisture chamber)

Wash grids with PBS

2nd antibody reaction

Room temp. 30min in gold-labeled 2nd antibody solution

Wash grids with DW

Electron staining

Uranium acetate 3min

Lead citrate 40sec

TEM observation