

Preparation of Crude Membrane Fraction

Homogenizing medium

0.25 M sorbitol	18.2 (g)
50 mM Tris/acetate (pH 7.5)*	0.5 M × 40 (mL)
1 mM EGTA*	0.1 M × 4 (mL)
1 mM MgCl ₂	0.1 M × 4 (mL)
(When you need to remove Mg ²⁺ , please add EDTA instead of Mg ²⁺)	
20 μM APMSF [#]	2 (mg)
1% (w/v) PVP	4 (g)
2 mM DTT	123 (mg)

Final volume	400 (mL)

[#] *p*-APMSF (Mr. 252.70): (*p*-amidinophenyl)methanesulfonyl fluoride hydrochloride (toxic)
APMSF should be added before experiment.

5% (w/v) Sucrose/Tris

Sucrose	5.0 g
20 mM Tris/acetate	0.5 M × 4 mL
1 mM EGTA	0.1 M × 1 mL
2 mM MgCl ₂	1 M × 0.2 mL
(When you need to remove Mg ²⁺ , please add EDTA instead of Mg ²⁺)	
2 mM DTT	31 mg

	100 mL

* 0.5 M Tris/acetate, pH 7.5 = 30.25 g / 500 mL
0.1 M EGTA = 19.0 g / 500 mL, adjust pH to 7.5 with KOH

Procedure

Plant Tissue (radish taproots or mung bean hypocotyls)

| ← homogenizing medium (chilled); [tissue : buffer = 1 : 1 or 1 : 2 (weight)]
| When you homogenize leaves, please add 5 to 10-fold volume of the buffer
| (10 g of leaves and 50 – 100 mL of buffer)

↓ homogenize with a Polytron

| 7,000 rpm × 10 min

Supernatant (Sup)

| 40,000 rpm × 25 min (RP45T, Beckman 45Ti)

Precipitate (Ppt)

| — suspend in 5% (w/v) Sucrose/Tris (or suitable buffer)

Crude membrane fractions

↓

Sucrose density gradient centrifugation (if you need)