

## Isolation of Plant Vacuolar Membranes

### Homogenizing medium

0.25 M sorbitol	18.2	36.4	54.7 (g)
50 mM Tris/acetate (pH 7.5)*	0.5 M soln. × 40	80	120 (mL)
1 mM EGTA*	0.1 M soln. × 4	8	12 (mL)
20 μM APMSF <sup>#</sup>	2	4	6 (mg)
1% (w/v) PVP	4	8	12 (g)
2 mM DTT	123	247	370 (mg)
<hr/>			
Final volume	400	800	1200 (mL)

<sup>#</sup> *p*-APMSF (Mr. 252.70): (*p*-amidinophenyl)methanesulfonyl fluoride hydrochloride (toxic)  
APMSF should be added before experiment.

### 0.5 M Sucrose/Tris

0.5 M sucrose	18 g
20 mM Tris/acetate	0.5 M × 4 mL
1 mM EGTA	0.1 M × 1 mL
2 mM MgCl <sub>2</sub>	1 M × 0.2 mL
2 mM DTT	31 mg

100 mL

### 0.25 M Sorbitol/Tris

0.25 M sorbitol	4.55 g
20 mM Tris/acetate	0.5 M × 4 mL
1 mM EGTA	0.1 M × 1 mL
2 mM MgCl <sub>2</sub>	1 M × 0.2 mL
2 mM DTT	31 mg

100 mL

\* 0.5 M Tris/acetate, pH 7.5 = 30.25 g / 500 mL  
0.1 M EGTA = 19.0 g / 500 mL, adjust pH to 7.5 with KOH

### Procedure

Plant Tissue (radish taproots or mung bean hypocotyls)
← homogenizing medium (chilled); [tissue : buffer = 1 : 1]
homogenize by a grater (Oroshigane in Japanese)
7,000 rpm × 10 min
<b>Supernatant (Sup)</b>
40,000 rpm × 25 min (RP45T, Beckman 45Ti)
<b>Precipitate (Ppt)</b>
← suspend in 0.5 M Sucrose/Tris
overlay with 0.25 M Sorbitol/Tris
40,000 rpm × 30 min (RP50T, Beckman 50.2Ti, slow brake)
<b>Interface between the two solutions</b>
← 0.25 M Sorbitol/Tris
40,000 rpm × 20 min
<b>Ppt</b>
← suspend in 0.25 M Sorbitol/Tris
(or 20 mM Tris-acetate/ 20% glycerol/1 mM DTT/
1 mM EGTA/1 mM MgCl <sub>2</sub> )
<b>Vacuolar membranes</b>