

Isolation of intact vacuoles from *Catharanthus* suspension-cultured cells

Preparation of protoplasts

In the case of about 24 g cells (Fresh weight from 160 ml suspension-cultured cells)

Remove the medium and wash the cells with the washing solution by aspiration

Solution A

Cellulase Y-C	1 g (1 %)
Pectolyase Y-23	50 mg (0.05 %)
P sor (0.5 M) pH 6.0	100 ml

↓

Incubate the cells in the solution A for 2.7 h at 31°C with shaking at 100 - 120 rpm

Purification of vacuoles

Move the released protoplasts to the 50 ml tube

Underlayered with Vc suc (0.4 M)

↓

Centrifugation (200 × g, 10 min)

↓

Remove the supernatant

↓

Add Vc suc

up to 15 ml and mix well

Form a gradient by overlaying

7 ml Vc sor (0.4 M) and

2 ml Vc GB (0.4 M)

↓

Centrifugation (200 × g, 1-2 min, and then 1600 × g, 8 min)

↓

Remove the solution over the purified protoplasts

Transfer the protoplasts to new tube




↓

Add equal volume of Vc med 0 and vortex

Incubate on the ice for 5 min

↓

Squeeze through a syringe 19 G (1.10 × 90 mm)



Add equal volume of VcGB and underlayered with Vc suc (0.4 M)

↓

Centrifugation (200 × g, 1-2 min, and then 1600 × g, 8 min)

↓

Remove the supernatant

↓

Divide the pellet into two glass tubes

↓

Add Vc suc (0.4 M) to half the tube

Form a gradient by overlaying

2 ml Vc suc : Vc sor = 1 : 2 ,

2 ml Vc sor, and

0.5 ml Vc GB

↓

Centrifugation (200 × g, 1-2 min, and then 1600 × g, 8 min)

↓

Vacuoles is obtained in the interphase (①)

between Vc GB and Vc sor

Collect solution (②)

After confirm vacuoles in the interphase (③),

Collect vacuoles

↓

Remove the layers over ④

↓

Squeeze through a syringe 19 G (1.10 × 90 mm) storonger than first time

↓

As described above, form a gradient

↓

Centrifugation (1000 rpm, 1-2 min, and then 3000 rpm, 8 min)

↓

Collect vacuoles




↓

In order to concentrate of vacuoles,

Add Vc suc (0.4 M) to

the solution of vacuoles, and overlay

Vc GB (0.4 M)






Centrifugation

(1000 rpm, 1-2 min, and then 3000 rpm, 8 min)



Collect vacuoles



Stock solutions

Protoplast 10 x	500 ml
Mes (213.25)	10.6625 g (100 mM)
CaCl ₂ · 2H ₂ O (147.02)	0.7351 g (10 mM)
pH 6.0 (Tris)	

Vc med 10 x	1000 ml
HEPES (238.3)	71.49 g (300 mM)
EGTA (380.4)	7.608 g (20 mM)
Potassium gluconate (234.2)	70.26 g (300 mM)
MgCl ₂ (203.3)	4.066 g (20 mM)
pH 7.2 with Tris	

Solutions

P sor (0.5 M) pH 6.0
10 mM Mes
1 mM CaCl ₂ · 2H ₂ O
0.5 M sorbitol

Washing solution

2 mM CaSO ₄ · 2H ₂ O
100 mM Sorbitol

Vc med 0

30 mM HEPES
2 mM EGTA
30 mM Potassium gluconate
2 mM MgCl ₂
pH 7.2 with Tris

Vc suc (0.4 M)	200 ml
Vc med 10 x	20 ml
2M sucrose	40 ml

Vc sor (0.4 M)	200 ml
Vc med 10 x	20 ml
2M sorbitol	40 ml

Vc GB (0.4 M)	200 ml
Vc med 10 x	20 ml
2M betaine monohydrate	40 ml

Vc suc : Vc sor = 1 : 2