


### Vanadate-Induced Nucleotide Trapping for Detection of ABC Proteins

1. Membrane proteins (20  $\mu\text{g}$ ) were incubated in a buffer mixture composed of 10  $\mu\text{M}$  8-azido-[ $\alpha$ - $^{32}\text{P}$ ]ATP (specific activity 740 GBq  $\text{mmol}^{-1}$ , ICN Biochemicals, Irvine, California), 2 mM ouabain, 0.1 mM EGTA, 3 mM  $\text{MgSO}_4$ , and 40 mM Tris-HCl (pH 7.5), 200  $\mu\text{M}$  orthovanadate in a total volume of 8  $\mu\text{l}$  for 10 min at 37°C.
2. The reactions were stopped by adding 400  $\mu\text{l}$  of ice-cold TGM buffer (40 mM Tris-HCl [pH 7.5], 0.1 mM EGTA, 1 mM  $\text{MgSO}_4$ ).
3. The membrane proteins were sedimented by centrifugation (20,000 $\times$ g, 10 min, 4°C) to separate unbound ATP.
4. After being washed with Tris-EGTA buffer, the pellet was irradiated with UV light at 254 nm (5.5  $\text{mW cm}^{-2}$ ) in 8  $\mu\text{l}$  of TGM buffer for 5 min on ice.
5. Samples were electrophoresed on a 7% SDS-polyacrylamide gel, and autoradiographed.
6. The radioactivity trapped by ABC proteins was analyzed with a radioimaging analyzer BAS 2000 (Fuji Film Co., Tokyo, Japan).



#### Reference

Kazuyoshi Terasaka, Nobukazu Shitan, Fumihiko Sato, Fumio Maniwa, Kazumitsu Ueda and Kazufumi Yazaki (2003)

Application of Vanadate-Induced Nucleotide Trapping to Plant Cells for Detection of ABC Proteins  
*Plant Cell Physiol.*, 44, 198-200.